

MENICARE SOFT



Indications

CLEANING • DEPROTEINIZATION • RINSING
DISINFECTING • STORAGE

How to use



1. Massage gently the lens for at least 20 seconds after depositing a few drops of MENICARE SOFT in the palm of your hand.
2. Rinse the lens thoroughly.
3. Soak for at least 4 hours.
4. Before wearing the lenses, rinse them with MENICARE SOFT.

Wide range of packaging



The MENICARE SOFT Care solution is a Class IIb medical device marked CE 0483, made in France, by Menicon Co., Ltd (regulatory manufacturer).
Read the instructions for use carefully.

New opening system*

Easy handling
Hygiene guaranteed



New ergonomic bottle*

Better grip
*To open the bottle, screw the cap clockwise as far as it will go. Note that this is only for 250 ml & 70 ml packaging.

MENICARE SOFT

Multi-purpose solution for all soft contact lenses



MP00090 2020-05 Photos: © all rights reserved - Printed in France by Menicon SAS - RCS Bobigny 722 038 668



Menicon Limited, Gatelodge Close, Round Spinney, Northampton NN3 8RJ
T +44 (0)1604 646 216 F +44 (0)1604 790 366 E enquiries@menicon.co.uk W www.menicon.co.uk

? > & * H # # ! \$ "

SAFETY COMFORT EFFICIENCY

For all soft contact lenses wearers




EFFECTIVE against germs, fungi & amibs
DEEP CLEANING of contact lenses


NON-CYTOTOXIC to corneal cells
BETTER comfort, hydration et lubrication

EFFICIENCY


- ✓ **Rapid removal** of germs and fungi
ISO 14729 Stand Alone Test¹: MENICARE SOFT meets the primary criteria for disinfecting solution in **20 minutes**.



- ✓ Over **99.9%** kill rate of *Acanthamoeba castellanii* (4 h)²
Against adhesion of *Anti-Acanthamoeba castellanii*³




- ✓ High **cleaning** capacity⁴
Macrolog HCO-60 is a **gentle and natural surfactant**



COMFORT

- ✓ **COMFORTEC**
A protective and hydrating film
Less protein deposits & Prevents protein denaturation⁵



Denatured protein deposited on the lens is correlated to ocular inflammation. Lysozyme is present in tears. In its original form, it protects eyes against germs. It denatures when it deposits on contact lenses. MENICARE SOFT prevents protein build-up and denaturation.⁵



SAFETY

MENICARE SOFT

- HARMONY BETWEEN COMFORTEC + PHMB
- NON-CYTOTOXIC¹
- PERFECT BALANCE

**High cleaning capacity
&
Broad spectrum disinfection**

**Respect for corneal cells
&
Good tolerance**

1.J.Santodomingo-Rubido, O.Mori, Satoshi Kawaminami, Cytotoxicity and antimicrobial activity of the new MeniCare Soft multi-purpose soft contact lens disinfecting solution. Ophthal. Physiol. Opt.2006 26(5): 476-472. 2.Lonnen J, Heaselgrave W, Nomachi M, et al. Disinfection efficacy and encystment rate of soft contact lens multipurpose solutions against Acanthamoeba. Eye & Contact Lens. 2010; 36:26-32. 3.Toshihiko Uno, MD, PhD, Yuichi Ohashi, MD, PhD, Miya Nomachi, MS, Masaki Imayasu, PhD Effects of Multipurpose Contact Lens Care Solutions on the Adhesion of Acanthamoeba to Silicone Hydrogel Contact Lenses. Cornea 2012;31:1170-1175. 4.Menicon R&D report n° C4685 / C5496 Cleaning efficacy on lipids. 5.Menicon R&D report n° C5233 and C5448 Protective efficacy against protein denaturation.