



Menicon Z Night

Clear vision during the day
without the need for glasses
or contact lenses

What is myopia?

Myopia, also known as near-sightedness, is an eye condition that is very widespread. If you have myopia, you can see things close to you clearly but objects that are far away are blurry.


Uncorrected myopia needs to be diagnosed and corrected to ensure you have optimal vision and can enjoy your life without visual restrictions.

Menicon Z Night lenses provide an effective and convenient alternative to spectacles, daily use contact lenses and refractive surgery.

Menicon Z Night lenses are worn at night and reshape your corneas whilst you sleep. The new corneal shape will correct your vision so when you remove the lenses every morning you will be free from needing lenses or spectacles during waking hours. This is a reversible procedure.



How Menicon Z Night works



Who is suitable for Menicon Z Night?

Menicon Z Night is for those with near-sightedness (myopia) up to -5.00D and astigmatism up to -4.00D.

Who can benefit from Menicon Z Night?

- Dissatisfied soft lens wearers
- People who want to be free from spectacles and daytime lens wear
- Sports enthusiasts who do not want the restrictions associated with lens wear during activities
- People looking for an alternative to refractive surgery

How can I get Menicon Z Night?

If you would like to know if you are a candidate to use Menicon Z Night, please contact your eye care professional. They will schedule an Initial Visit to check your vision, prescription, eye measurements and eye health to ensure you are suitable and will order your trial lenses.

From your side, all you need is some time for the fitting process and the motivation!

**Feel the freedom of not needing any
visual correction during the day**

Maintenance

Your Menicon Z Night lenses will be used every night for one year. Then, they will need to be replaced.

To keep your lenses in the best condition: Menicon recommends MeniCare Pure (or Plus) for daily lens cleaning, rinsing, and disinfection, and intensive cleaning with Menicon Progent regularly.

For further details, please see the instructions for use for these products.





Menicon Ltd.

T: +44 (0)1604 646216

E: enquiries@menicon.co.uk

www.menicon.co.uk

Menicon Z Night is manufactured by Menicon B.V. (CE0344),
MeniCare Plus, MeniCare Pure and Progent are manufactured by
Menicon Pharma SAS (CE0483) and are medical devices.

MZNBCBUKENP.01